





DHRS9 Antibody

Product Code	CSB-PA006867DA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BPW9
Immunogen	Recombinant Human Dehydrogenase/reductase SDR family member 9 protein (23-319AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dehydrogenase/reductase SDR family member 9 (EC 1.1.-.-) (3-alpha hydroxysteroid dehydrogenase) (3-alpha-HSD) (NADP-dependent retinol dehydrogenase/reductase) (RDH-E2) (RDHL) (Retinol dehydrogenase 15) (Short chain dehydrogenase/reductase family 9C member 4) (Short-chain dehydrogenase/reductase retSDR8) (Tracheobronchial epithelial cell-specific retinol dehydrogenase) (RDH-TBE), DHRS9, RDH15 SDR9C4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DHRS9

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA006867DA01HU at dilution of 1:100



Immunofluorescent analysis of A549 cells using CSB-PA006867DA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)