




rpiA Antibody

Product Code	CSB-PA540419DA01ENT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A7Z0
Immunogen	Recombinant Escherichia coli Ribose-5-phosphate isomerase A protein (1-219AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribose-5-phosphate isomerase A (EC 5.3.1.6) (Phosphoriboisomerase A) (PRI), rpiA, ygfC
Immunogen Species	Escherichia coli
Research Area	Signal Transduction
Target Names	rpiA

Image



Western blot
 All lanes: rpiA antibody at 2µg/ml + DH5a whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 23 kDa
 Observed band size: 23 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.