



# esxB Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA363658YC01MVZ
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P9WNK4
<b>Immunogen</b>	Recombinant Mycobacterium tuberculosis ESAT-6-like protein EsxB protein (2-100AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mycobacterium tuberculosis
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	ESAT-6-like protein EsxB (10 kDa culture filtrate antigen CFP-10) (Secreted antigenic protein MTSA-10), esxB, cfp10 lhp mtsA10
<b>Immunogen Species</b>	Mycobacterium tuberculosis (strain CDC 1551 / Oshkosh)
<b>Research Area</b>	Others
<b>Target Names</b>	esxB