



ATP4B Antibody

Product Code	CSB-PA002343EA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51164
Immunogen	Recombinant Human Potassium-transporting ATPase subunit beta protein (58-291AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Potassium-transporting ATPase subunit beta (Gastric H(+)/K(+) ATPase subunit beta) (Proton pump beta chain), ATP4B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ATP4B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.