




# Tpsb2 Antibody

<b>Product Code</b>	CSB-PA024128LA01MO
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P21845
<b>Immunogen</b>	Recombinant Mouse Tryptase beta-2 protein (32-276AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mouse, Rat
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:500-1:2000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Tryptase beta-2 (Tryptase-2) (EC 3.4.21.59) (Mast cell protease 6) (mMCP-6), Tpsb2, Mcpt6
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Research Area</b>	Immunology
<b>Target Names</b>	Tpsb2

## Image



Western blot  
 All lanes: Tpsb2 antibody at 2µg/ml  
 Lane 1: Mouse heart tissue  
 Lane 2: Rat heart tissue  
 Secondary  
 Goat polyclonal to rabbit IgG at 1/10000 dilution  
 Predicted band size: 31, 26 kDa  
 Observed band size: 31 kDa