



IL12B Antibody, Biotin conjugated

Product Code	CSB-PA011587YD01SH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68220
Immunogen	Recombinant Sheep Interleukin-12 subunit beta protein (23-327AA)
Raised In	Rabbit
Species Reactivity	Sheep
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-12 subunit beta (IL-12B) (Cytotoxic lymphocyte maturation factor 40 kDa subunit) (CLMF p40) (IL-12 subunit p40), IL12B
Immunogen Species	Ovis aries (Sheep)
Research Area	Immunology
Target Names	IL12B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.