



# Recombinant Dog Somatotropin (GH1)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-EP009407DO   |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.  |
| <b>Uniprot No.</b>       | P33711   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Canis lupus familiaris (Dog) (Canis familiaris)  |
| <b>Purity</b>            | ≥85% (SDS-PAGE)  |
| <b>Sequence</b>          | FPAM PLSSLFANAV LRAQHLHQLA ADTYKEFERA YIPEGQRYSI<br>QNAQA AFCFS ETIPAPTGKD EAQQRSDVEL LRFSLLIQS WLGPVQFLSR<br>VFTNSLVFGT SDRVYEKLD LEEGIQALMR ELEDGSPRAG QILKQTYDKF<br>DTNLRSDDAL LKNYGLLS CF KKDLHKAETY LRVMKCRRFV ESSCAF   |
| <b>Source</b>            | E.coli   |
| <b>Target Names</b>      | GH1  |
| <b>Protein Names</b>     | Recommended name: Somatotropin Alternative name(s): Growth hormone   |
| <b>Expression Region</b> | 27-216   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | Full Length of Mature Protein  |
| <b>Target Details</b>    | This protein is a member of the somatotropin/prolactin family of hormones which play an important role in growth control. The gene, along with four other related genes, is located at the growth hormone locus on chromosome 17 where they are interspersed in the same transcriptional orientation; an arrangement which is thought to have evolved by a series of gene duplications. The five genes share a remarkably high degree of sequence identity. Alternative splicing generates additional isoforms of each of the five growth hormones, leading to further diversity and potential for specialization. This particular family member is expressed in the pituitary but not in placental tissue as is the case for the other four genes in the growth hormone locus. Mutations in or deletions of the gene lead to growth hormone deficiency and short stature. |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |