



# Recombinant Chicken Protein BTG1 (BTG1)

<b>Product Code</b>	CSB-EP002863CH-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P34743
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Gallus gallus (Chicken)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MHPALYTRAS MIREIAAAVA FISKFLRTKG LMNERQLQTF SQSLQELLAE HYKHHWFPEK PCKGSGYRCI RINHKMDPLI GQAAQRIGLS SQELFQLLPS ELTLWVDPYE VSYRIGEDGS ICVLYEAAPA GGSQNNTNMQ MVDSRISCKE ELLLGRTSPS KSYNMMTVSG
<b>Source</b>	E.coli
<b>Target Names</b>	BTG1
<b>Protein Names</b>	Recommended name: Protein BTG1 Alternative name(s): B-cell translocation gene 1 protein
<b>Expression Region</b>	1-170
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	This gene is a member of an anti-proliferative gene family that regulates cell growth and differentiation. Expression of this gene is highest in the G0/G1 phases of the cell cycle and downregulated when cells progressed through G1. The encoded protein interacts with several nuclear receptors, and functions as a coactivator of cell differentiation. This locus has been shown to be involved in a t(8;12)(q24;q22) chromosomal translocation in a case of B-cell chronic lymphocytic leukemia.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.