



Recombinant Human Glypican-1 (GPC1)

Product Code	CSB-BP009703HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P35052
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	DPASKSR SCGEVRQIYG AKGFSLSDVP QAEISGEHLR ICPQGYTCCT SEMEENLANR SHAELETALR DSSRVLQAML ATQLRSFDDH FQHLLNDSER TLQATFPGAF GELYTQNARA FRDLYSELRL YYRGANLHLE ETLAEFWARL LERLFKQLHP QLLLPDDYLD CLGKQAEALR PFGEAPREL RLRATRAFVAA RSFVQGLGVA SDVVRKVAQV PLGPECSRVA MCLVYCAHCL GVPGARPCPD YCRNVLKGCL ANQADLDAEW RNLLDSMVLITDKFWGTSGV ESIVIGSVHTW LAEAINALQD NRDTLTAKVI QGCGNPKVNP QGPGPEEKRR RGKLAPRERP PSGTLEKLVS EAKAQLRDVQ DFWISLPGTL CSEKMALSTA SDDRCWNGMA RGRYLPEVMG DGLANQINNP EVEVDITKPD MTIRQQIMQL KIMTNRLRSA YNGNDVDFQD ASDDGSGSGS GDGCLDDLCS RKVSRKSSSS RTPLTHALPG LSEQEGQKTS
Source	Baculovirus
Target Names	GPC1
Protein Names	Recommended name: Glypican-1 Cleaved into the following chain: 1. Secreted glypican-1
Expression Region	24-530
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	Cell surface heparan sulfate proteoglycans are composed of a membrane-associated protein core substituted with a variable number of heparan sulfate chains. Members of the glypican-related integral membrane proteoglycan family (GRIPS) contain a core protein anchored to the cytoplasmic membrane via a glycosyl phosphatidylinositol linkage. These proteins may play a role in the control of cell division and growth regulation.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.