




Recombinant Human CMRF35-like molecule 1 (CD300LF), partial

Product Code	CSB-YP819470HU
Relevance	Acts as an inhibitory receptor for myeloid cells and mast cells. Inhibits osteoclast formation.
Abbreviation	Recombinant Human CD300LF protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TDQ1
Alias	CD300 antigen-like family member FImmune receptor expressed on myeloid cells 1 ;IREM-1Immunoglobulin superfamily member 13 ;IgSF13NK inhibitory receptor; CD300f
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	TQITGPTTVNGLERGLSTVQCVYRSGWETYLKWWCRGAIWRDCKILVKTS SGSEQEVKRDVRVSIKDNQKNRTFTVTMEDLMKTDADTYWCGIEKTGNDLGVT VQVTIDPAPVTQEETSSSPTLTGHHLNDRHKLKLS
Research Area	Immunology
Source	Yeast
Target Names	CD300LF
Protein Names	Recommended name: CMRF35-like molecule 1 Short name= CLM-1 Alternative name(s): CD300 antigen-like family member F Immune receptor expressed on myeloid cells 1 Short name= IREM-1 Immunoglobulin superfamily member 13 Short
Expression Region	20-156aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	17.5kDa
Protein Length	Extracellular Domain
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.