




# Recombinant Mouse Serine protease inhibitor A3N (Serpina3n)

<b>Product Code</b>	CSB-YP835578MO
<b>Relevance</b>	The single human alpha1-antichymotrypsin gene (SERPINA3) is represented by a cluster of 14 individual murine paralogs.
<b>Abbreviation</b>	Recombinant Mouse Serpina3n protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q91WP6
<b>Storage Buffer</b>	Tris-based buffer,50% glycerol
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	FPDGTLGMDAAVQEDHDNGTQLDSLTLASINTDFAFSLYKELVLKNPDKNIVFS PLSISAALAVMSLGAKGNTLEEILEGLKFNLTETSEADIHQQFGHLLQRLNQPK DQVQISTGSALFIEKRQQILTEFQEAKACTLYQAEAFADFDQQRQAKKLINDYV RKQTQGMIKELVSDLDKRTLMLVNVYIYFKAKWKVPFDPLDTFKSEFYAGKRR PVIVPMMSMEDLTTPYFRDEELSCTVVELKYTGNASALFILPDQGRMQQVEAS LQPETLRKWKNSLKPRMIDELHLPKFSISTDYSLEDVLSKLGIREVVFSTQADLSA ITGTKDLRVSQVVHKAVLDVAETGTEAAAATGVKFPMSAKLYPLTVYFNRPFLL IMIFDTETEIAPFIAKIANPK
<b>Research Area</b>	Others
<b>Source</b>	Yeast
<b>Target Names</b>	Serpina3n
<b>Protein Names</b>	Spi2
<b>Expression Region</b>	21-418aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	46.8 kDa
<b>Protein Length</b>	Full Length of Mature Protein
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.