



L-lactate dehydrogenase Antibody

Product Code	CSB-PA632935LA01PLS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q27743
Immunogen	Recombinant Plasmodium falciparum L-lactate dehydrogenase protein (238-316AA)
Raised In	Rabbit
Species Reactivity	Plasmodium falciparum
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	L-lactate dehydrogenase (EC 1.1.1.27) (LDH-P)
Immunogen Species	Plasmodium falciparum
Research Area	Others