




# CALCA Monoclonal Antibody

<b>Product Code</b>	CSB-MA1213441A0m
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P01258
<b>Immunogen</b>	Recombinant Human PCT
<b>Raised In</b>	mouse
<b>Species Reactivity</b>	Human
<b>Specificity</b>	specific for Human PCT denatured and native forms
<b>Tested Applications</b>	ELISA,IHC;Recommended dilution:IHC:1:50-1:500
<b>Form</b>	liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%,protein G purified
<b>Clonality</b>	monoclonal
<b>Alias</b>	PCT
<b>Product Type</b>	Monoclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	CALCA
<b>Clone No.</b>	4G1H7

**Image**


Immunohistochemical of paraffin-embedded human thyroid tissues using CSB-MA1213441A0m at dilution of 1:200

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.