




PGC Monoclonal Antibody

Product Code	CSB-MA112991A0m
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20142
Immunogen	Recombinant Human PGC
Raised In	mouse
Species Reactivity	Human
Specificity	specific for PGC native and denatured forms
Tested Applications	ELISA,IHC;Recommended dilution:IHC:1:50-1:500
Relevance	Hydrolyzes a variety of proteins.
Form	liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%,protein G purified
Isotype	IgG1
Clonality	monoclonal
Alias	Pepsinogen C ,PGC
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PGC
Clone No.	1G1H3

Image



Immunohistochemical of paraffin-embedded human gastric cancer tissue using CSB-MA112991A0m at dilution of 1:200

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.