



SERPINB1 Antibody, HRP conjugated

Product Code	CSB-PA021065LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30740
Immunogen	Recombinant Human Leukocyte elastase inhibitor protein (1-379AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leukocyte elastase inhibitor (LEI) (Monocyte/neutrophil elastase inhibitor) (EI) (M/NEI) (Peptidase inhibitor 2) (PI-2) (Serpine B1), SERPINB1, ELANH2 MNEI PI2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SERPINB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.