




ZSCAN20 Antibody

Product Code	CSB-PA027136LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17040
Immunogen	Recombinant Human Zinc finger and SCAN domain-containing protein 20 protein (1-433AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger and SCAN domain-containing protein 20 (Zinc finger protein 31) (Zinc finger protein 360) (Zinc finger protein KOX29), ZSCAN20, KOX29 ZNF31 ZNF360
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZSCAN20

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA027136LA01HU at dilution of 1:100