



EML2 Antibody, HRP conjugated

Product Code	CSB-PA007643LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95834
Immunogen	Recombinant Human Echinoderm microtubule-associated protein-like 2 protein (1-418AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Echinoderm microtubule-associated protein-like 2 (EMAP-2) (HuEMAP-2), EML2, EMAP2 EMAPL2
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	EML2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.