



DENR Antibody

Product Code	CSB-PA006724HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43583
Immunogen	Recombinant Human Density-regulated protein (2-198AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Density-regulated protein (DRP) (Protein DRP1) (Smooth muscle cell-associated protein 3) (SMAP-3), DENR, DRP1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	DENR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.