



SENp8 Antibody

Product Code	CSB-PA822251HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96LD8
Immunogen	Recombinant Human Sentrin-specific protease 8 protein (1-212AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sentrin-specific protease 8 (EC 3.4.22.-) (Deneddylase-1) (NEDD8-specific protease 1) (Protease, cysteine 2) (Sentrin/SUMO-specific protease SENP8), SENP8, DEN1 NEDP1 PRSC2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SENp8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.