



LCK Antibody, HRP conjugated

Product Code	CSB-PA012806LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06239
Immunogen	Recombinant Human Tyrosine-protein kinase Lck protein (1-539AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosine-protein kinase Lck (EC 2.7.10.2) (Leukocyte C-terminal Src kinase) (LSK) (Lymphocyte cell-specific protein-tyrosine kinase) (Protein YT16) (Proto-oncogene Lck) (T cell-specific protein-tyrosine kinase) (p56-LCK), LCK
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	LCK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.