




# RAB11A Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA019153HA01HU   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                    |
| <b>Uniprot No.</b>         | P62491   |
| <b>Immunogen</b>           | Recombinant Human Ras-related protein Rab-11A protein (2-213AA)                  |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Human  |
| <b>Tested Applications</b> | ELISA, IHC; Recommended dilution: IHC:1:20-1:200                                 |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | Non-conjugated   |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| <b>Purification Method</b> | >95%, Protein G purified   |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Alias</b>               | Ras-related protein Rab-11A (Rab-11) (YL8), RAB11A, RAB11                        |
| <b>Immunogen Species</b>   | Homo sapiens (Human)   |
| <b>Research Area</b>       | Signal Transduction  |
| <b>Target Names</b>        | RAB11A   |

## Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA019153HA01HU at dilution of 1:100