



NTMT1 Antibody

Product Code	CSB-PA880111HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BV86
Immunogen	Recombinant Human N-terminal Xaa-Pro-Lys N-methyltransferase 1 protein (1-223AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	N-terminal Xaa-Pro-Lys N-methyltransferase 1 (EC 2.1.1.244) (Alpha N-terminal protein methyltransferase 1A) (Methyltransferase-like protein 11A) (N-terminal RCC1 methyltransferase) (X-Pro-Lys N-terminal protein methyltransferase 1A) (NTM1A) [Cleaved into: N-terminal Xaa-Pro-Lys N-methyltransferase 1, N-terminally processed], NTMT1, C9orf32 METTL11A NRMT NRMT1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NTMT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.