



IDUA Antibody, FITC conjugated

Product Code	CSB-PA011000LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35475
Immunogen	Recombinant Human Alpha-L-iduronidase protein (28-653AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-L-iduronidase (EC 3.2.1.76), IDUA
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	IDUA