



CBLIF Antibody, Biotin conjugated

Product Code	CSB-PA009417LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27352
Immunogen	Recombinant Human Cobalamin binding intrinsic factor protein (19-417AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cobalamin binding intrinsic factor (Gastric intrinsic factor) (Intrinsic factor) (IF) (INF), CBLIF, GIF, IFMH
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CBLIF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.