





ELANE Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA007587LA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P08246 |
| Immunogen | Recombinant Human Neutrophil elastase protein (30-267AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Neutrophil elastase (EC 3.4.21.37) (Bone marrow serine protease) (Elastase-2) (Human leukocyte elastase) (HLE) (Medullasin) (PMN elastase), ELANE, ELA2 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | ELANE |

Image



Immunohistochemistry of paraffin-embedded human spleen tissue using CSB-PA007587LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA007587LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.