



ELANE Antibody, HRP conjugated

Product Code	CSB-PA007587LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08246
Immunogen	Recombinant Human Neutrophil elastase protein (30-267AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	Modifies the functions of natural killer cells, monocytes and granulocytes. Inhibits C5a-dependent neutrophil enzyme release and chemotaxis.
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neutrophil elastase (EC 3.4.21.37) (Bone marrow serine protease) (Elastase-2) (Human leukocyte elastase) (HLE) (Medullasin) (PMN elastase), ELANE, ELA2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	ELANE