





KYAT1 Antibody

Product Code	CSB-PA613699LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16773
Immunogen	Recombinant Human Kynurenine--oxoglutarate transaminase 1 protein (1-372AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kynurenine--oxoglutarate transaminase 1 (EC 2.6.1.7) (Cysteine-S-conjugate beta-lyase) (EC 4.4.1.13) (Glutamine transaminase K) (GTK) (Glutamine--phenylpyruvate transaminase) (EC 2.6.1.64) (Kynurenine aminotransferase 1) (Kynurenine aminotransferase I) (KATI) (Kynurenine--oxoglutarate transaminase I), KYAT1, CCBL1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KYAT1

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA613699LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA613699LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.