



Genome polyprotein Antibody, HRP conjugated

Product Code	CSB-PA18549B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17763
Immunogen	Recombinant Dengue virus type 1 Genome polyprotein protein
Raised In	Rabbit
Species Reactivity	Dengue virus type 1
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Genome polyprotein [Cleaved into: Capsid protein C (Capsid protein) (Core protein); Protein prM (Precursor membrane protein); Peptide pr (Peptide precursor); Small envelope protein M (Matrix protein); Envelope protein E; Non-structural protein 1 (NS1); Non-structural protein 2A (NS2A); Serine protease subunit NS2B (Flavivirin protease NS2B regulatory subunit) (Non-structural protein 2B); Serine protease NS3 (EC 3.4.21.91) (EC 3.6.1.15) (EC 3.6.4.13) (Flavivirin protease NS3 catalytic subunit) (Non-structural protein 3); Non-structural protein 4A (NS4A); Peptide 2k; Non-structural protein 4B (NS4B); RNA-directed RNA polymerase NS5 (EC 2.1.1.56) (EC 2.1.1.57) (EC 2.7.7.48) (Non-structural protein 5)]
Immunogen Species	Dengue virus type 1
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.