



AMH Antibody, FITC conjugated

Product Code	CSB-PA12639C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03971
Immunogen	Recombinant Human Mullerian-inhibiting factor protein (26-560AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mullerian-inhibiting factor (Anti-Mullerian hormone) (AMH) (Mullerian-inhibiting substance) (MIS), AMH, MIF
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AMH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.