



KCNJ10 Antibody, FITC conjugated

Product Code	CSB-PA012048LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78508
Immunogen	Recombinant Human ATP-sensitive inward rectifier potassium channel 10 protein (165-379AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-sensitive inward rectifier potassium channel 10 (ATP-dependent inwardly rectifying potassium channel Kir4.1) (Inward rectifier K(+) channel Kir1.2) (Potassium channel, inwardly rectifying subfamily J member 10), KCNJ10
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCNJ10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.