



# AKT2 Antibody, Biotin conjugated

Product Code	CSB-PA15919D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31751
Immunogen	Recombinant Human RAC-beta serine/threonine-protein kinase protein (2-481AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RAC-beta serine/threonine-protein kinase (EC 2.7.11.1) (Protein kinase Akt-2) (Protein kinase B beta) (PKB beta) (RAC-PK-beta), AKT2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AKT2