



NPC2 Antibody, FITC conjugated

Product Code	CSB-PA634688LC01DO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q28895
Immunogen	Recombinant Dog NPC intracellular cholesterol transporter 2 protein (22-149AA)
Raised In	Rabbit
Species Reactivity	Dog
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NPC intracellular cholesterol transporter 2 (Epididymal secretory protein E1) (cE1) (Niemann Pick type C2 protein homolog), NPC2
Immunogen Species	Dog
Research Area	Others
Target Names	NPC2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.