



Cd68 Antibody

Product Code	CSB-PA004951LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31996
Immunogen	Recombinant Mouse Macrosialin protein (21-291AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Macrosialin (CD antigen CD68), Cd68
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Cd68