



Cd68 Antibody, FITC conjugated

Product Code	CSB-PA004951LC01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31996
Immunogen	Recombinant Mouse Macrosialin protein (21-291AA)
Raised In	Rabbit
Species Reactivity	Mouse
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Macrosialin (CD antigen CD68), Cd68
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Cd68
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.