




Pdia2 Antibody

Product Code	CSB-PA017719HA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	D3Z6P0
Immunogen	Recombinant Mouse Protein disulfide-isomerase A2 protein (21-527AA)
Raised In	Rabbit
Species Reactivity	Mouse, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein disulfide-isomerase A2 (EC 5.3.4.1) (PDIp), Pdia2, Pdip
Immunogen Species	Mus musculus (Mouse)
Research Area	Tags & Cell Markers
Target Names	Pdia2

Image



Western blot
 All lanes: Pdia2 antibody at 10µg/ml + Hela whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 59, 58 kDa
 Observed band size: 59 kDa

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.