



Ttr Antibody

Product Code	CSB-PA025270HA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07309
Immunogen	Recombinant Mouse Transthyretin protein (23-147AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transthyretin (Prealbumin), Ttr
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Ttr
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.