



Gstp1 Antibody, Biotin conjugated

Product Code	CSB-PA009989HD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19157
Immunogen	Recombinant Mouse Glutathione S-transferase P 1 protein (2-210AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutathione S-transferase P 1 (Gst P1) (EC 2.5.1.18) (GST YF-YF) (GST class-pi) (GST-piB) (Preadipocyte growth factor), Gstp1, Gstpib
Immunogen Species	Mus musculus (Mouse)
Research Area	Cardiovascular
Target Names	Gstp1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.