



Fc μ r Antibody, HRP conjugated

Product Code	CSB-PA007972LB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A1KXC4
Immunogen	Recombinant Mouse Fas apoptotic inhibitory molecule 3 protein (18-262AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fas apoptotic inhibitory molecule 3 (IgM Fc fragment receptor) (Regulator of Fas-induced apoptosis Toso), Fc μ r, Faim3 Toso
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Fc μ r
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.