



# Fc $\mu$ r Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA007972LC01MO
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	A1KXC4
<b>Immunogen</b>	Recombinant Mouse Fas apoptotic inhibitory molecule 3 protein (18-262AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Fas apoptotic inhibitory molecule 3 (IgM Fc fragment receptor) (Regulator of Fas-induced apoptosis Toso), Fc $\mu$ r, Faim3 Toso
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Research Area</b>	Others
<b>Target Names</b>	Fc $\mu$ r
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.