



rgpA Antibody, Biotin conjugated

Product Code	CSB-PA338957LD01PQP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28784
Immunogen	Recombinant Porphyromonas gingivalis Gingipain R1 protein (228-720AA)
Raised In	Rabbit
Species Reactivity	Porphyromonas gingivalis
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gingipain R1 (EC 3.4.22.37) (Arg-gingipain) (Gingipain 1) (RGP-1), rgpA, rgp1
Immunogen Species	Porphyromonas gingivalis
Research Area	Others
Target Names	rgpA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.