



CXCL10 Antibody

Product Code	CSB-PA822646LA01MOW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8MIZ1
Immunogen	Recombinant Macaca mulatta C-X-C motif chemokine 10 protein (22-98AA)
Raised In	Rabbit
Species Reactivity	Macaca mulatta
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-X-C motif chemokine 10 (10 kDa interferon gamma-induced protein) (Gamma-IP10) (IP-10) (Small-inducible cytokine B10), CXCL10, SCYB10
Immunogen Species	Macaca mulatta (Rhesus macaque)
Research Area	Others
Target Names	CXCL10