



CXCL10 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA822646LA01MOW |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8MIZ1 |
| Immunogen | Recombinant Macaca mulatta C-X-C motif chemokine 10 protein (22-98AA) |
| Raised In | Rabbit |
| Species Reactivity | Macaca mulatta |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | C-X-C motif chemokine 10 (10 kDa interferon gamma-induced protein) (Gamma-IP10) (IP-10) (Small-inducible cytokine B10), CXCL10, SCYB10 |
| Immunogen Species | Macaca mulatta (Rhesus macaque) |
| Research Area | Others |
| Target Names | CXCL10 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |