



S Antibody

Product Code	CSB-PA336163EA01HIY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36334
Immunogen	Recombinant Human coronavirus OC43 Spike glycoprotein protein (15-344AA)
Raised In	Rabbit
Species Reactivity	Human coronavirus OC43
Tested Applications	ELISA
Relevance	S1 attaches the virion to the cell membrane by interacting with sialic acid-containing cell receptors, initiating the infection. S2 is a class I viral fusion protein. Under the current model, the protein has at least 3 conformational states: pre-fusion native state, pre-hairpin intermediate state, and post-fusion hairpin state. During viral and target cell membrane fusion, the coiled coil regions (heptad repeats) assume a trimer-of-hairpins structure, positioning the fusion peptide in close proximity to the C-terminal region of the ectodomain. The formation of this structure appears to drive apposition and subsequent fusion of viral and target cell membranes.
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Spike glycoprotein (S glycoprotein) (E2) (Peplomer protein) [Cleaved into: Spike protein S1 (90B); Spike protein S2 (90A)], S
Immunogen Species	Human coronavirus OC43
Research Area	Others
Target Names	S