



# inhA Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA363781HC01MVZ
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P9WGR1
<b>Immunogen</b>	Recombinant Mycobacterium tuberculosis Enoyl-[acyl-carrier-protein] reductase [NADH] protein (1-269AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mycobacterium tuberculosis
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Enoyl-[acyl-carrier-protein] reductase [NADH] (ENR) (Enoyl-ACP reductase) (EC 1.3.1.9) (FAS-II enoyl-ACP reductase) (NADH-dependent 2-trans-enoyl-ACP reductase), inhA
<b>Immunogen Species</b>	Mycobacterium tuberculosis (strain ATCC 25618 / H37Rv)
<b>Research Area</b>	Others
<b>Target Names</b>	inhA
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.