



inhA Antibody, Biotin conjugated

Product Code	CSB-PA363781HD01MVZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P9WGR1
Immunogen	Recombinant Mycobacterium tuberculosis Enoyl-[acyl-carrier-protein] reductase [NADH] protein (1-269AA)
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Enoyl-[acyl-carrier-protein] reductase [NADH] (ENR) (Enoyl-ACP reductase) (EC 1.3.1.9) (FAS-II enoyl-ACP reductase) (NADH-dependent 2-trans-enoyl-ACP reductase), inhA
Immunogen Species	Mycobacterium tuberculosis (strain ATCC 25618 / H37Rv)
Research Area	Others
Target Names	inhA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.