




Serpib2 Antibody

Product Code	CSB-PA021070HA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29524
Immunogen	Recombinant Rat Plasminogen activator inhibitor 2 type A protein (1-416AA)
Raised In	Rabbit
Species Reactivity	Rat, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Plasminogen activator inhibitor 2 type A (PAI-2) (PAI2A) (Serpib B2), Serpib2, Pai2 Planh2
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Others
Target Names	Serpib2

Image



Western blot
 All lanes: Serpib2 antibody at 2µg/ml
 Lane 1: MCF-7 whole cell lysate
 Lane 2: 293T whole cell lysate
 Lane 3: Rat liver tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 48 kDa
 Observed band size: 48 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.