



Serpinb2 Antibody

Product Code	CSB-PA021070HA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29524
Immunogen	Recombinant Rat Plasminogen activator inhibitor 2 type A protein (1-416AA)
Raised In	Rabbit
Species Reactivity	Rat, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Plasminogen activator inhibitor 2 type A (PAI-2) (PAI2A) (Serpin B2), Serpinb2, Pai2 Planh2
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Others
Target Names	Serpinb2

Image



Western blot

All lanes: Serpinb2 antibody at 2μg/ml

Lane 1: MCF-7 whole cell lysate

Lane 2: 293T whole cell lysate

Lane 3: Rat liver tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 48 kDa

Observed band size: 48 kDa