



# Genome polyprotein Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA362073LD01HQD
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P07210
<b>Immunogen</b>	Recombinant Human rhinovirus A serotype 89 Genome polyprotein protein (575-866AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human rhinovirus A serotype 89
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Genome polyprotein [Cleaved into: P3; Protein 3AB; P2; P1; Capsid protein VP0 (VP4-VP2); Capsid protein VP4 (P1A) (Virion protein 4); Capsid protein VP2 (P1B) (Virion protein 2); Capsid protein VP3 (P1C) (Virion protein 3); Capsid protein VP1 (P1D) (Virion protein 1); Protease 2A (P2A) (EC 3.4.22.29) (Picornain 2A) (Protein 2A); Protein 2B (P2B); Protein 2C (P2C) (EC 3.6.1.15); Protein 3A (P3A); Viral protein genome-linked (VPg) (Protein 3B) (P3B); Protein 3CD (EC 3.4.22.28); Protease 3C (P3C) (EC 3.4.22.28); RNA-directed RNA polymerase (RdRp) (EC 2.7.7.48) (3D polymerase) (3Dpol) (Protein 3D) (3D)]
<b>Immunogen Species</b>	Human rhinovirus A serotype 89
<b>Research Area</b>	Others
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.