



Single-stranded DNA-binding protein Antibody, Biotin conjugated

Product Code	CSB-PA366021HD01EEB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03696
Immunogen	Recombinant Enterobacteria phage T7 Single-stranded DNA-binding protein gp2.5 protein (1-232AA)
Raised In	Rabbit
Species Reactivity	Enterobacteria phage T7
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterobacteria phage T7
Research Area	Others
Target Names	2.5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.