



Single-stranded DNA-binding protein Antibody, Biotin conjugated

Product Code	CSB-PA366021HD01EEB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03696
Immunogen	Recombinant Enterobacteria phage T7 Single-stranded DNA-binding protein gp2.5 protein (1-232AA)
Raised In	Rabbit
Species Reactivity	Enterobacteria phage T7
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterobacteria phage T7
Research Area	Others
Target Names	2.5